



Technical Data Sheet

Release date 01/06/2024

Product Name: Exo Elixir Liquorice		Product code: COS_EEL	
Product Description: Spray-dried extract of Liquorice (<i>Glycyrrhiza glabra</i> L., dried root)			
INCI Name / CAS number / EC number of Ingredients: GLYCYRRHIZA GLABRA EXTRACT / CAS number: 84775-66-6 / EC number: 283-895-2 MALTODEXTRIN / CAS number: 9050-36-6/ EC number: 232-940-4			
Main Specifications			
Composition:	Spray-dried extract of: [Liquorice (<i>Glycyrrhiza glabra</i> L., dried root); Maltodextrin]		
Origin of vegetables:	Organic agriculture		
Country of origin:	Italy		
pH value at 25 ± 0.2°C:	5.6 (measured diluted 1:10 in MilliQ® water)		
Recommended period of use:	Within 24 months from production if stored under optimal conditions.		
Mode of use:	See attached document "Product use guidelines" for details.		
Guidance on Storage:	Store the product at room temperature (20 - 25°C), tightly closed, away from heat sources, not exposed to light and humidity.		
Storage and preservation requirements			
Storage temperature:	At room temperature (20 - 25°C) in spray-dried form (See "Guidance on Storage" section for details)		
Transport temperature:	At room temperature (20 - 25°C)		
Organoleptic Properties			
Color:	Ivory white		
Appearance:	Spray-dried powder		
Consistence:	Spray-dried powder		
Chemical-Physical Characteristics		Microbiological Characteristics	
Humidity	< 5 %MC Test method: Rapporti ISTISAN 1996/34 pag 7 Met. A	Total combined Yeast-Mould Count	< 10 CFU/g Test method: ISO 16212
Ash	< LOQ Test method: Rapporti ISTISAN 1996/34 pag 77	Total Aerobic Mesophilic Bacteria Count	< 10 CFU/g Test method: ISO 21149



Heavy metals	Cadmium (Cd): < 0.06 ± 0.03 mg/Kg Mercury (Hg): < 0.21 ± 0.08 mg/kg Lead (Pb): < 0.25 ± 0.10 mg/kg Test method: UNI EN 13805:2014 + UNI EN 15763:2010	Candida albicans	Absent Test method: ISO 18416
Polycyclic Aromatic Hydrocarbons	< LOQ Test method: UNI EN 16619:2015	Escherichia coli	Absent Test method: ISO 21150
Multiresidual analysis of pesticides	< LOQ Test method: UNI EN 15662:2018	Pseudomonas aeruginosa	Absent Test method: ISO 22717
		Staphylococcus aureus	Absent Test method: ISO 22718
		Enterobacteriacee	< LOQ Test method: UNI EN ISO 21528-2:2017/EC 1:2018
		Salmonella spp	Absent Test method: UNI EN ISO 6579-1:2020
		Aflatoxins B1, B2, G1, G2	< LOQ Test method: MI D 0015 rev.01 2023
Compliant with EC Regulation 2073/2005 and successive amendments and supplements			

Statement on GMOs (EC Reg. 1829/2003 and EC Reg. 1830/2003)
The product does not contain GMOs and is not consisting of substances containing GMOs
Delivery time
45 days from the order
Packing and transport requirements
Primary packaging: glass containers
Secondary packaging: isothermal polystyrene containers
Tertiary packaging: carton containers
Abbreviations
LOQ: Limit of Quantification



exolab
I T A L I A

The Laboratory Manager

Dr. Mariantonia Logozzi

Mariantonia Logozzi